



NEW!



ClearNote HD

Transform your laptop into a powerful magnifier! The ClearNote HD connects to your computer and provides you with magnified, high contrast live images of your surroundings, such as presentation handouts and white boards.

The ClearNote HD is easy to use. It is touchscreen compatible with your laptop or tablet, so you can simply pinch to zoom, swipe to move and touch to select. No need to remember complicated key combinations.

Get sharp and clear images, without compromising color quality. With the automatic contrast enhancement blurry print is enhanced while viewing the images and text in color.

An ultra-flexible camera head rotates both vertically and horizontally, providing a crisp and clear live image from any distance. To quickly change your focus from a document to a white board, previously selected settings of different camera positions are automatically recalled, making reading quicker and easier.

Conveniently designed to draw power from a computer there is no need for an extra power supply.





Optelec ClearNote HD

Key benefits



- Small and compact with easy set-up, perfect for on-the-go
- Combines near and distance viewing in one design
- Automatic contrast enhancement for sharp, clear images
- Save and store images on computer for later viewing
- Touchscreen support
- High Definition camera for superior images
- Continuous magnification: 1.3X to 40X (on a 15" monitor)
- Powered by USB
- Super-light, weighs only 3.1 lbs
- Connects with computer, laptop or Windows tablet (Windows 8 or higher) via USB 3.0

Available through your Authorized Optelec Dealer:

Optelec U.S. Inc.

Tel: 800.826.4200 | www.Optelec.com

**3030 Enterprise Court, Suite C
Vista, CA 92081**

Optelec
Life is worth enjoying



magnification



speech



daily living aids

Optelec U.S. Inc.

Tel: 800.826.4200 - Fax: 800.368.4111 - info@optelec.com - www.optelec.com

Optelec has offices in:

THE NETHERLANDS – GERMANY – BELGIUM – UNITED KINGDOM – USA – CANADA

Simple and effective solutions for the blind and visually impaired